

**ELECTRONIC BLOOD** PRESSURE MONITOR

Shenzhen Pango Medical Electronics Co., Ltd Main Site: Building 2, No. 25 Fenghuang Road, Industrial Zone, Xikeng First Village, Henggang Street, Longgang District, Shenzhen, 518115 Guangdong, P. R. China. Additional site I: 2-4 Floor ,No.5 Shanzhuang Rd., Xikeng Village, Henggang Town, Longgang District, Shenzhen, 518115 Guangdong, P. R. China. Tel:+86-755-33825988 Fax:+86-755-33825989

EC REP Lotus NL B.V. Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. Tel: +31644168999

Date: 2023-04-13 Rev:A/2

**TABLE OF CONTENTS** INTRODUCTION..... NOTES ON SAFETY......2 PRECAUTIONS BEFORE USE......8 FEATURES OF THE PRODUCT......9 PARTS IDENTIFICATION......10 INSERT OR REPLACE BATTERIES......11 TIME AND VOICE ON/OFF OF SYSTEM SETUP......12 UNIT CONVERSION mmHg/kPa DISPLAY......13 WHO BLOOD PRESSURE CLASSIFICATION DISPLAY......14 HOW TO MEASURE BLOOD PRESSURE......16 CLEAN AND MAINTENANCE......17 CALIBRATION METHOD......29

INTRODUCTION The Monitor uses the oscillometric method of blood pressure measurement.

Measurement Automatic Electronic Blood Pressure Monitor is intended for use by medical professionals or at home to monitor and display diastolic. systolic blood pressure and pulse rate, with an air wrist cuff buckled around one's wrist according to the instructions in the "ATTACHING THE WRIST CUFF." The expected life of the product is 5 years. The product complies with the electromagnetic compatibility requirement of IEC 60601-1-2 and safety standards of IEC 60601-1 and performance of

IEC 80601-2-30 as specified in Regulation (EU)2017/745. **NOTES ON SAFETY** 

\* The warning signs and sample icons shown here are listed for your safe and correct use of the unit, so as to prevent injuries or damages to the device. \* The icons and meanings are as follow.

Examples of signs The  $\otimes$  icon indicates prohibitions (what you should not do). Matters involving actual prohibitions are indicated by text or pictures in or near . The left icon refers to "general prohibition".

C E 2862

The • icon indicates something that is compulsory (what must always Patient must follow doctor's instruction and should not perform self-judgment and self-treatment by the measuring result, Matters involving actual compulsory actions are indicated by text or Self-diagnosis of measured results and treatment are dangerous, pictures in or near . The left icon refers to "general compulsion". The device should not be used to judge illness, first aid and continuously monitor measuring. The So icon indicates something can't be disassembled or "Don' disassemble" Matters involving actual compulsory actions are indicated by text or pictures

is MR-unsafe:

∕!\ Caution Type BF Applied part **IP Classification: IP20 Consult instruction** The following Please refer to the for use symbol indicates instructions for use that the device

in or near 

. The left icon refers to general prohibition.

Indicates a medical device that needs to be protected from moisture.

5

Do not smoke

Reduce salt and fat intake

3. Why measure blood pressure at home?

Maintain proper weight

be observed)

MD Indicates

Contact its local authorities to determine the proper method

of disposal of potentially bio hazardous parts and accessories.

Requests from Manufacturer Make sure there is no connection tubing kinking before start measuring to

Exercise regularly

Have regular physical checkups

medical device

This device can not be used for Patient transport and surgical care .It can be used in household or fixed places only. Please press "on/off" button to stop work when you feel uncomfortable with the wrist, or if the air is inflating abnormally without stop. This device should not be used by children under 18 years old or people who cannot express their will, otherwise it will cause harm. Do not use the unit for purpose other than measuring blood pressure.

May cause accident or trouble. Please do not use mobile phone around the device. Please do not use the device around the magnetic field. The device is prohibited from being used during movement.

Do not use the equipment in outdoor or shower rooms. Do not disassemble, repair, or remodel the main unit or the wrist cuff of the blood pressure monitor. Will cause the unit to function erroneously

-Not servicing and maintenance while the ME EQUIPMENT is in use.

-Stop using the equipment immediately, if it is in contact with water.

the maintenance instructions of manual

hypertensive individuals, variations are

Normally, the blood pressure rises while

even more pronounced.

-The user can maintain the product, the maintenance method is described in

-The PATIENT is an intended OPERATOR

6

Typical fluctuation within a day

(Measured every five minutes)

2

Į.

avoid any injury to patient. For any patient, do not measure more than 3 times continuously, it should be at least above 5 minutes of interval rest between any two measurements, otherwise will cause extravasated blood. Do not measure your blood pressure over 6 times each day. Do not apply the cuff over a wound as this can cause further injury. Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury. Observe the air pressure value from the LCD display. When measuring, it could not exceed 280 mmHg, otherwise Please press "on/off" button to stop Do not use force to bend the wrist cuff or the air tube. Do not knock or drop the main unit. Always use the specified accessories in the manual, the use of other parts not approved by the manufacturer may cause faults or injuries For service information, parts list etc., please contact the dealer.

**ABOUT BLOOD PRESSURE** 1. What is blood pressure? Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing

in the morning while one is still at rest and before eating. 2. What is hypertension and how is it controlled? Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.

To prevent hypertension or keep it under control:

Blood pressure measured at a clinic or doctor's office may cause apprehension concerned by the results of one and produce an elevated reading, 25 to 30 mmHg higher than that measured at home, Home measurement reduces the effects of outside influences on blood measurement. pressure readings, supplements the doctor's readings and provides a more Take measurements at the same time every accurate, complete blood pressure history. day using the procedure described in this manual, and know your normal blood pressure. 4. WHO blood pressure classification Many readings give a more comprehensive Standards for assessment of high blood blood pressure history. pressure, without regard to age, have

Reference Material: Journal of Hypertension been established by the World Health 1999, Vol 17 No.2 Organization (WHO), and shown in Grade 3 hypertension (severe) 2 110 Grade 2 hypertension (moderate) 5. Blood pressure variations 월 100 Grade 1 hypertension (mild 95 An individual's blood pressure varies 90 greatly on a daily and seasonal basis. High-normal

It may vary by 30 to 50 mmHg due to various conditions during the day. In

85 Normai Optimal 120 130 140 150 160 170 180 Systolic blood pressure mr 2. For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be fluctuation in blood pressure values measured at the upper arm versus at the wrist.

at work or play and falls to its lowest levels during sleep. So, do not be overly Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.

**PRECAUTIONS BEFORE USE** 

1. If you are taking medication, consult with your doctor to determine the most

appropriate time to measure your blood pressure. NEVER change a prescribed

medication without first consulting with your doctor.

3. WHO blood pressure classification display.

Systolic Blood pressure Diastolic Blood pressure

• The batteries may leak and cause a malfunction.

months or more), remove the batteries.

testing monitor performance and may have a shorter life.

nomal boot unit values are shown as blood pressure.

Also select memory unit value changes.

mmHc

150

8

4. Easy to use, Press a button to automatically measure, record the measurement

values and measurement time. 5. Automatically turns off (within 1 minute) to save power. PARTS IDENTIFICATION SYMBOLS ON DISPLAY LCD Display

devices or turn them off. 4. Before using, should wash your hands. 5. Do not measure on the arm which simultaneously used monitoring ME Equipment, otherwise it could cause loss of function. 6. Consult your doctor if the unexpected readings are obtained, also please refer to "Trouble shooting" of the manual. 7. The reading is probably a little lower than measured in the hospital due to the steady mood at home.

3. Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields.

To prevent such interference, use the monitor at a sufficient distance from such

8.Cuff pressure range 0-299mmHg **FEATURES OF THE PRODUCT** 1. Memory can store 90 measurements. 2. Large and clear LCD display.

**INSERT OR REPLACE BATTERIES** 

2. Insert new batteries into the battery compartment as shown, taking care that the

1. Remove the battery cover.

polarities(+) and (-)are correct.

model does not have this function.

Wrist Cuf

WHO blood press SET Button ON/OFF Button Accessory Manual 10 · Batteries, which have fluid on surface or be modified, can not be inserted into the products. Battery short circuit must be prevented. Battery life varies with the ambient temperature and may be shorten at low

• Use the specified batteries only. The batteries provided with the device are for

Used batteries may leak and damage the main unit. Pleases observe the following

\* If you are not going to use the unit for a long period of time (approximately three

## Disposal of empty battery to the authorized collecting party

3. Close the battery cover, Use only LR03, AAA batteries.

subject to the regulation of each individual territory. CAUTION

• Insert the batteries as shown in the battery compartment. If not, the device will not work.

ullet When lacksquare (LOW BATTERY mark) blinks in the display, replace all batteries with new

ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction. (LOW BATTERY mark) does not appear when the batteries run out. • Please ensure to distinguish positive polar "+" and negative polar "-" of batteries when replacing batteries. 11

4. Press "MEM" key to adjust the month. Following the same steps to adjust date/hour/

minute/Voice (on/off) until setting completed ("  $\square \Pi$ " is the On, "  $\square P$ " is the Off) Non-talking

TIME AND VOICE ON/OFF OF SYSTEM SETUP 1. Press "SET" key to Time display. 2. In the off state, Press and hold "SET" key until the year number displays and flashes on LCD to enter setting mode. 3. Press "MEM" key to adjust the year, then press "SET" key

\* Replace worn batteries with their polarities in the correct direction.

again to save your setting and enter the month setting mode.

The units will be chosen by the above shows mmHg/kPa after decontrol, After the

year month date hour minute Voice UNIT CONVERSION mmHg/kPa DISPLAY The goods have mm Hg(mmHg), kPa (kPa) two kinds of blood pressure display units(mmHg factory to express).

Press "ON / OFF" button for 10 seconds to display unit switching interface, then

ATTACHING THE WRIST CUFF

1) Wrap the wrist cuff around your wrist about (1-2)cm above your hand as shown

press "MEM" key to select mmHg / KPa, press "ON / OFF" button to exit.

Reference material: journal of

**"**-\ondo

ÖΠ

WHO BLOOD PRESSURE CLASSIFICATION DISPLAY Grade 3 hypertension (severe) Diastolic blood pressure

hypertension 1999. vol 17 No.2

measurement, and display "0" or last measurement record

3. Start measurement, the cuff in the strap will automatically inflate.

The mark(♥) will flash on LCD. When complete, the results will be displayed.

**HOW TO MEASURE BLOOD PRESSURE** 1. Fasten the wrist cuff according to the instructions in "ATTACHING THE WRIST CUFF." 2. Press the "ON/OFF" button. All icons appear two seconds on DISPLAY, then switch to

Grade 2 hypertension (moderate)

14

16

18

20

Grade 1 hypertension (mild)

High-normal

Normal Optimal

in the figure at the right. 2) Fasten the wrist cuff tightly by using the Velcro Strip. For proper measurements, fasten the wrist cuff tightly and measure on a bare wrist. 2. How to take proper measurements

15

1. Fastening the wrist cuff

13

## • Relax for about 5 to 10 minutes before measurement. • Raise your hand so that the wrist cuff is at the same

level as your heart. • Remain still and keep quiet during measurement. • Do not measure left after physical exercise or a bath. Measure your blood pressure at about the same time every day.

**READ MEMORY** 

Measuring Method

Measuring Range:

Indication

Accuracy:

Power supply:

Operating condition:

Storage condition:

Dimensions:

Classification

Wrist circumference

Weight:

Memory:

For best accuracy in blood pressure measurement: • Sit comfortably at a table. Rest your wrist on the table.

Press " MEM " button to inquire memory average values " RU9 "Average Value Display: the latest 3 groups of memory average values (Memory values are displayed regardless of period). Press "MEM" button, a memory reading out the latest measurements, "MEM" for the buttons(UP), "SET" button for the memory (DOWN) Power Measurement closure or after the end of the state .can press the "MEM" button read out the latest measurement of memory. **DELETE MEMORY** The state read out the memory press the (memory) button five seconds, the LCD display "П□" has been to delete all memory

∞0: 68 **CLEAN AND MAINTENANCE** 1. Keep this device in the case provided with the device when you do not use it. 2.Do not fold the arm cuff too tightly.

**SPECIFICATIONS** 

Digital LCD display

90 Memories

Type BF

(13.5~19.5)cm \* Specifications may be changed without notice in the event of improvement being made.

Oscillometric Measurement

Pressure:(30~280)mmHg Pulse:(40~199)Beat/min

2x1.5V Batteries(LR03 or AAA)

+5°C~+40°C. 15%RH~93%RH

-20°C~+55°C. 0%RH~93%RH

Approx: 76(W)X68(H)X32(D)mm

Approx: 120g, excluding batteries

Static Pressure:  $\pm 3$ mmHg Pulse:  $\pm 5\%$ 

Atmospheric pressure: 70kPa~106kPa

Atmospheric pressure:50kPa~106kPa

use alkaline battery, measure above 200 times.

ПΟ

3. Clean the monitor with a soft dry cloth. Do not use any cleaning solution. 4.Do not submerge the device or any components in water 5. Store the device and the components in a clean and safe location. 6. The clean steps for the cuff is provided as following. \* Completely wipe the inner side (the side that contacts skin) of the cuff with a soft cloth lightly moistened with 75% Ethyl alcohol 3 times. Replace the soft cloth after each wipe. \*Then air dry the cuff.

\* Remove the batteries if the unit will not be used for three months or longer. Always replace all the batteries with new ones at the same time • This product is designed for use over an extended period of time; however, it is generally recommended that it be inspected and calibrated every two years to ensure proper function and performance. \* See the Calibration Method for more details.

Storage conditions: -20°C~+55°C. 0%RH~93%RH

**ERROR DISPLAY** 

Nothing is displayed When you push the

POWER button or

**Emissions** 

RF emissions

RF emissions

CISPR 11

CISPR 11

Battery icon flash

that it is used in such an environment.

Compliance

Group 1

Class B

ISM bandsa

\* Do not submerge the device or any of the components in water. Do not subject the monitor to extreme hot or cold temperatures,

\* Store the device and the components in a clean, safe location.

**CAUTION** 

humidity or direct sunlight.

1. Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT. 2.Degree or protection against electric shock: TYPE BF APPLIED PART. 3. Mode of operation: CONTINUOUS OPERATION. 4. Equipment not suitable for category AP&APG equipment use in presence. **STATEMENT** the system might not meet its performance specifications if stored or used outside the temperature and humidity as mentioned below: Operating conditions: +5°C~+40°C. 15%RH~93%RH 70kPa~106kPa

TROUBLE SHOOTING

**HOW TO CORRECT** 

Replace new batteries

Insert battery in the correct

Insert batteries

Electromagnetic environment-guidance

The Model PG-800A25 Series Electronic Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions

The Model PG-800A25 Series Electronic Blood

including cables, than the recommended

separation distance calculated from the

equation applicable to the frequency of

Recommended separation distance

24

the transmitter.

are very low and are not likely to cause any interference in nearby electronic equipment

Pressure Monitor is used in home and it's

polarities

If you have trouble in using the unit please check the following points first.

POSSIBLE CAUSE

Battery worn out

placed wrongly

**Appendix 1 Guidance and Manufacturer Declaration Tables** 


Guidance and manufacturer's declaration – electromagnetic emissions The Model PG-800A25 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A25 Series Electronic Blood Pressure Monitor should assure

No battery installation

The polarities of batteries

Do not subject the monitor to strong shocks, such as dropping the unit on the floor

E1:can't normally Increase pressure	Check your wrist cuff if any air leakage	Replace wrist cuff with new one
E3 inflate pressure too high	Pressure value of more than 299mmHg	Re-measurement or send bac dealer for re-calibrate pressure
E2E4:have shaking while measurement	Hand or body shaking while measurement	keeping static and correct gesture to measure again
Battery icon on	Battery low power	Replace battery and measure agai
The systolic pressure Value or diastolic Pressure value too high	1.The wrist cuff was held lower than your heart	keeping correct position and gesture to measure again
	2.The wrist cuff was not attached properly	
	3. You moved your body or spoke during measurement	
The systolic pressure Value or diastolic Pressure value too low	1.The wrist cuff was held higher than your heart	
	2.you moved your body or Spoke during measurement	



Model PG-800A25 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment IEC 60601 Compliance Immunity test Electromagnetic environmenttest level quidance level Electrostatic ±8 kV contact ±8 kV contact Floors should be wood, concrete ±2 kV, ±4 kV, ±2 kV, ±4 kV, or ceramic tile. If floors are covered discharge ±8 kV. ±15KV  $\pm 8$  kV,  $\pm 15$  KV with synthetic material, the relative (ESD)IEC air humidity should be at least 30 %. 61000-4-2 Power 30 A/m, 50/60Hz 30 A/m, 50/60Hz Power frequency magnetic fields should be at levels frequency characteristic of a typical (50/60 Hz) location in a typical commercial magnetic field or hospital environment. IEC 61000-4-8 NOTE  $U_{\tau}$  is the a.c. mains voltage prior to application of the test level 23 Radiated RF 10 V/m  $d = \sqrt{\frac{3.5}{E_1}} \sqrt{P} \quad 80 \text{MHz to } 800 \text{MHz}$ 10 V/m IEC 61000-4-3 80 MHz to 2.7 GHz

 $d = \left| \frac{7}{E_1} \right| \sqrt{P}$  800MHz to 2.7GHz

where P is the maximum output power

and d is the recommended separation

according to the transmitter manufacturer

Field strengths from fixed RF transmitters,

as determined by an electromagnetic site

compliance level in each frequency range

Interference may occur in the vicinity

survey, a should be less than the

rating of the transmitter in watts (W)

distance in metres(m).

Guidance and manufacturer's declaration - electromagnetic immunity

The Model PG-800A25 Series Electronic Blood Pressure Monitor is intended for use

in the electromagnetic environment specified below. The customer or the user of the

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7.3 MHz. 10.1 MHz to 10.15 MHz. 14 MHz to 14.2 MHz. 18.07 MHz to 18.17 MHZ, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50.0 MHz to 54.0 MHz b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges. 26

of equipment marked with the following symbol: (((•))) 25 c Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-800A25 Series Electronic Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Model PG-800A25 Series Electronic Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-800A25 Series Electronic Blood Pressure Monitor. d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. Recommended separation distances between portable and mobile RF communications equipment and the Model PG-800A25 **Series Electronic Blood Pressure Monitor** 

For transmitters rated at a maximum output power not listed above the

recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum

output power rating of the transmitter in watts (W) according to the transmitter

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency

27

29

range applies.

The Model PG-800A25 Series Electronic Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model PG-800A25 Series Electronic Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model PG-800A25 Series Electronic Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment. Rated maximum Separation distance according to frequency of transmitter output of transmitter 150 kHz to 80 MHz | 80 MHz to 800 MHz | 800 MHz to 2.7 GHz  $d = \left[\frac{3.5}{V_1}\right] \sqrt{P}$  $d = \left[\frac{7}{E_1}\right] \sqrt{P}$ W

> 0.12 0.38

1.2

10 3.8 3.8 7.3 100 12 12 23 28 4. External input 50mmHg and 200mmHg standard static air pressure, and observe the air pressure value displayed at the position of the LCD systolic pressure (SYS) and the value of the digital pressure gauge should be in the range of +/-3mmHg.

0.38

1.2

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

1. Press and hold the "ON/OFF, MEM" button at the same time, load the battery, enter the static air pressure calibration mode after the LCD screen is fully displayed, and then release the button. 2. Press ON/OFF to close the internal air valve. 3. Connect the external standard barometric interface and the digital barometer

**CALIBRATION METHOD** 

interface to the cuff interface.

interference environment to ensure basic safety and basic performance unchanged. 2.In the event of any serious event related to this product, such as serious adverse

0.01

0.1

1

Essential performance: Limits of the error of the manometer, ±3mmHg.Reproducibility

102

Clinical benefits: Accurate measurement of SBP and DBP, clinical performance meets

30

0.73

2.3

<u>∕!</u>\ Caution 1. ME devices can be used in exposed environments, including electromagnetic event, significant alteration of the product resulting in change of intended use, etc., it will be reported to the manufacturer and the competent authorities of the user and/or the member states where the patient is located. of the blood pressure determination, ±3mmHg. the requirements of ISO 81060-2:2018.